

ActiveCare™ Bed

**Safely enter
and exit the bed
from a full chair
position with the
push of ONE button.**



Since the inception of the hospital bed - one of the biggest challenges for caregivers and patients is **safely** getting in and out of the bed. Today, Med-Mizer has innovated the **ActiveCare** Bed to solve this challenge.

- **One Button SafeTurn Exit** — exclusive feature that fully automates the bed to turn into a chair and allow the patient to safely enter and exit the bed. Use the built in high/low as a stand assist to offer additional assistance and safety.
- **Full Chair Position** — full chair position allows the patient the ability to sit in bed for added comfort. Watch TV, eat a meal, and socialize with family and friends while remaining comfortably in bed.
- **RollBack in Place™** — head of bed can articulate to almost 90 degrees with the assistance of rollback in place to alleviate shearing/friction in the lower back.
- **Home-Like Features** — enjoy wood aesthetics and hidden casters, optional night light, optional USB charging to provide your patients the ultimate in "home-like" comfort and amenities
- **Advanced Features** — optional built in scale, out of bed exit alarm/light, ability to sync into nurse call, digital HOB, trendelenburg to create the most advanced "residential looking" hospital bed

ActiveCare® Bed Deluxe Model SPECIFICATIONS

Deck Height Low	13.25"
Deck Height High	28.5"
Sleeping Surface	35"w, 42"w x 80"l
Maximum Safe Working Capacity	600lbs


LIMITED WARRANTY

15 years on bed frame mechanics
2 years on drive system


1 year on wood, casters &
handsets
Lifetime warranty on all welds




ENDBOARD OPTIONS



Obsidian Oak




Gray Oak




Natural Walnut

MATTRESS OPTIONS



Deluxe - 7" Multi Layered foam with middle die cut for superior pressure redistribution



Elite - 7" Copper gel infused memory foam for superior immersion and heat dissipation

SIDERAIL OPTIONS



Mobility Bar



Telescoping

ACTIVECARE PORTFOLIO



Standard - Powered Height with OneButtonSafeTurn



Fixed - Fixed Height with OneButtonSafeTurn

PROUDLY MADE IN THE
U.S.A.

(877) 867-7365
MED-MIZER.COM

MED MIZER®

is at the forefront of innovating medical technology for post acute care. We are focused on partnering with you to help address healthcare challenges in the following key areas—




FALLS
PREVENTION



BARIATRIC
CARE



MOBILITY



REHABILITATION



WOUND CARE

ActiveCare™ Series



Fixed Height

- ✓ One Button Safe Turn
- ✓ Fully Adjustable Head and Knee Sections

Optional Features:

- ✓ 3" or 5" casters
- ✓ Battery BackUp
- ✓ Mobility Bar
- ✓ Up/Down Side Rail



Standard

- ✓ One Button Safe Turn
- ✓ Built in Stand Assist (High/Low Height Adjustment)
- ✓ Fully Adjustable Head and Knee Sections

Optional Features:

- ✓ Trend/Reverse Trend Pendant
- ✓ Trapeze
- ✓ Battery BackUp
- ✓ Mobility Bar
- ✓ Up/Down Side Rail



Deluxe

- ✓ One Button Safe Turn
- ✓ Built in Stand Assist (High/Low Height Adjustment)
- ✓ Fully Adjustable Head and Knee Sections
- ✓ Premium Aesthetics with hidden casters
- ✓ Telescoping Rails
- ✓ Battery BackUp and Nightlight

Optional Features:

- ✓ Trend/Reverse Trend Pendant
- ✓ Touchscreen with Lockouts
- ✓ Mobility Bar