

ATTO

Mobility never looked so good



See the world – and be seen – with ATTO, the smartest high-performance mobility scooter on the market. While the compact folded unit splits in two for easier transport, ATTO comes together to form a robust full-sized vehicle with the best balance and use of space.

ATTO's beauty is in the details, with standard features like a built-in USB charging port, deck storage and adjustable seat heights. It can also be customized with accessories like hidden arm rests and batteries with various capacities. A head-turning piece of equipment for those who literally would not ride in anything else, ATTO is the mobility scooter that refuses to look like one.

See for yourself at www.movinglife.com

DIMENSIONS

Trolley mode: H 28" 71 cm
L 15.3" 38.9 cm
W 16.5" 41.9 cm

Drive mode: H 35" 89 cm
L 47" 120 cm
W 22" 56 cm

WEIGHT

Front half: 26.4 lbs 12 kg
Rear half: 35.7 lbs 16 kg
Assembled: 62.2 lbs 28.2kg (w/out battery)
Battery: 4.1 lbs 1.85 kg

WEIGHT CAPACITY

Maximum user weight:
220 lbs 100 kg

BATTERY

48 Volts lithium-ion, more than 10 miles of continuous driving, only 4 hours on average to fully charge, charge on-board and off-board

SPEED

Most countries, 4 mph/ 6.4 kilometer per hour (ATTO can go up to 6 mph)

Separates into 2 parts
in seconds

Easy to wheel
as a trolley

ATTO
Mobility never looked so good



movinglife™

Compact to
full size in **10** sec