

# BEDS

CLINICAL FUNCTIONALITY, HOME-LIKE STYLE



## KEY FACTS

- › 9.8" Low Height
- › 13" Transfer Position
- › 32" High Height
- › Sleep Surface: 39" x 79", External Dimensions: 44" x 81"
- › Extra-Wide Frame Upgrade: Sleep Surface: 48" x 79" External Dimensions: 53" x 81"
- › Space-saving, vertical side rail
- › Four-wheel central brake
- › Integrated bed extension
- › Four-section reclining surface
- › Automatic comfort seating mechanism
- › Safe working load of (500 lbs)

## FEATURES




Adjustable comfort

The adjustable mattress deck height provides support while getting in and out of bed, customizable to the height and preference of the user.




Non-institutional, elegant design

Various finish options compliment any room decor, promoting comfort and tranquility for users and caregivers.



The Comfort Deck surface flexes to provide ultimate comfort to the user.

Each section of the frame can be easily removed for thorough cleaning.




Illuminating hand controller

The bed frame is adjusted with a hand controller that is conveniently located at the side of the bed. Simple images on the pendant indicate which buttons control various sections of the bed. Buttons illuminate in the dark for easy visibility.

## ONE-TOUCH POSITIONING

The hand controller features various pre-set positioning options including Chair, Auto Contour, Mobility Assistance, Trendelenburg, Reverse Trendelenburg and more, giving the user security and comfort at the touch of a button.



- + Mobility Assistance
- + Independence
- + Comfort

## WOODGRAIN/UPHOLSTERY SELECTION




Slate Grey

SonderCare beds are available in two finishes to match your modern or classic design style.



Slate Grey



Espresso

## MATTRESS SELECTION

The SonderCare pressure relieving mattress is designed to work hand-in-hand with your bed system to provide optimal comfort and enjoyment.

- + Cool Memory Foam
- + Durable Foundation Layer
- + Two Cover Styles
  - 2-Way Stretch, fluid-proof \*not currently available on the Extra-Wide frame
  - Dream Quilt Bamboo Cover
- + Dimensions
  - 39" x 79" (Twin XL), 48" x 79" (Extra-Wide Upgrade)
  - 6" high (Stretch Cover) or 7" high (Quilt Cover)



Contact us directly for more information  
or to schedule a visit to our showroom  
[www.sondercare.com](http://www.sondercare.com) // [info@sondercare.com](mailto:info@sondercare.com) // 833.649.7772

**SonderCare**  
Healthy Living Adjustable Beds 

# MATTRESSES

DESIGNED TO MEET YOUR NEEDS FOR THE SLEEP OF YOUR DREAMS

## Comfort & Dream Memory Foam

The SonderCare pressure relieving mattress is designed to work hand-in-hand with your bed system to provide optimal comfort and enjoyment. There are two memory foam mattress options to choose from, our Comfort or Dream option.



*SonderCare Dream  
Mattress shown above*



*SonderCare Comfort  
Mattress shown above*



*SonderCare Dream Mattress shown above on  
the Aura™ Premium Bed*



### KEY FACTS

- › Cool Memory Foam
- › Durable Foundation Layer

#### Two Cover Styles

- › 2-Way, Stretch, fluid-proof (comes standard on the Comfort Mattress)
- › Padded Quilt Cover (comes standard on the Dream Mattress)

#### Dimensions

- › 39" x 80" (Twin XL)
- › Comfort Mattress: 6" depth  
Dream Mattress: 7" depth

