

MADE IN  
**USA**  
OF FOREIGN & DOMESTIC COMPONENTS

# VivaLIFT!®

POWER RECLINERS

By Pride® Mobility Products Corp.

# THE MOST COMFORTABLE RECLINER IN THE WORLD!

## A TRUE-INFINITE 4-MOTOR RECLINER

WITH  
COUNTLESS  
POSITION  
COMBINATIONS



Shown in Durasoft Oat Fabric



# ESCAPE



Enjoy endless relaxation with the Escape Collection from **VivaLift®** Power Recliners. Power lumbar offers personalized comfort and a full-width power headrest lifts your head, neck and shoulders. A convenient remote with USB charger lets you power your smart device, and the industry-first Lithium battery backup ensures your lift chair functions when the power goes out. With an array of gorgeous fabrics, it's an easy decision to make this fabulous collection a part of your home.


**Pride Mobility®** (US) 800-800-8586 | (CAN) 888-570-1113  
pridemobility.com | prideliftchairs.com |    

# ESCAPE



COLLECTION	ESCAPE	ESCAPE
MODEL	PLR-990iM	PLR-990iL
POSITION	True-Infinite Position Zero Gravity Trendelenburg Chaise Lounger	True-Infinite Position Zero Gravity Trendelenburg Chaise Lounger
SIZE	Medium	Large
WEIGHT CAPACITY	375 lbs.	375 lbs.
SEAT-TO-FLOOR	20"	20"
SEAT DEPTH	20"	22"
SEAT WIDTH	20"	23"
TOP OF BACK TO SEAT	26"	29"
FULLY PADDED CHAISE	Yes	Yes
FOOTREST EXTENSION	Standard	Standard
SEATING TYPE	Pocketed Coil Springs	Pocketed Coil Springs
CHAIR HEIGHT	44.5"	47"
DISTANCE FROM WALL	18"	21"
OVERALL WIDTH	33"	36.5"
CHAIR WEIGHT	148 lbs.	150 lbs.
WARRANTY	7-Yr. Pro-rated	7-Yr. Pro-rated

## FUNCTIONS AND FEATURES:


-  POWER HEADREST
-  POWER LUMBAR
-  USB REMOTE
-  LITHIUM BATTERY BACKUP
-  FOOTREST EXTENSION
-  POWER BACKREST

## POSITION GUIDE:


-  RELAXED
-  RECLINED
-  UPWARD LIFT
-  TRUE INFINITE



## STANDARD FABRIC: DURASOFT



OAT



TIMBER

 OPTIONAL FABRIC:  
**ULTRALEATHER™** 



BUFF



FUDGE



Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



= Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

**Pride Mobility®** (US) 800-800-8586 | (CAN) 888-570-1113

pridemobility.com | prideliftchairs.com |    