

VivaLIFT![®]
POWER RECLINERS
By Pride Mobility Products[®] Corp.

THE MOST
COMFORTABLE
RECLINER
IN THE WORLD!

THE 4-MOTOR RECLINER

WITH
COUNTLESS
POSITION
COMBINATIONS



Shown in Astro
Brown Fabric

TRANQUIL









There is no greater joy than reading a good book in a comfortable, cozy recliner. The Tranquil Collection by **VivaLift!**[®] Power Recliners combines quiet comfort with gorgeous fabrics, for a rich and upscale look. Get the support you need for your neck, shoulders and lower back as you dive into your favorite story. The Tranquil Collection exudes style and class with contrast stitching. It's the perfect addition to your home.





TRANQUIL

COLLECTION	TRANQUIL	TRANQUIL	TRANQUIL	TRANQUIL
MODEL	PLR-935PW	PLR-935S	PLR-935M	PLR-935LT
POSITION	Infinite Lay Flat Relaxed Reclined Upward Lift	Infinite Lay Flat Relaxed Reclined Upward Lift	Infinite Lay Flat Relaxed Reclined Upward Lift	Infinite Lay Flat Relaxed Reclined Upward Lift
SIZE	Petite Wide	Small	Medium	Large/Tall
WEIGHT CAPACITY	375 lbs.	375 lbs.	375 lbs.	375 lbs.
SEAT-TO-FLOOR	18"	18"	19.5"	21.5"
SEAT DEPTH	18.5"	18.5"	20.5"	22.5"
SEAT WIDTH	24"	20"	22"	23.5"
TOP OF BACK TO SEAT	24.5"	25"	27"	30"
FULLY PADDED CHAISE	Yes	Yes	Yes	Yes
FOOTREST EXTENSION	Standard	Standard	Standard	Standard
SEATING TYPE	Total Comfort	Total Comfort	Total Comfort	Total Comfort
CHAIR HEIGHT	38"	40"	42"	47"
DISTANCE FROM WALL	12.5"	13"	16"	16"
OVERALL WIDTH	36"	33.75"	37.5"	38"
CHAIR WEIGHT	143 lbs.	137 lbs.	145 lbs.	155 lbs.
WARRANTY	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated

FUNCTIONS AND FEATURES:

-  POWER HEADREST
-  POWER LUMBAR
-  USB REMOTE
-  LITHIUM BATTERY BACKUP
-  FOOTREST EXTENSION
-  POWER BACKREST

POSITION GUIDE:

-  RELAXED
-  RECLINED
-  UPWARD LIFT
-  INFINITE LAY FLAT

STANDARD FABRIC: ASTRO



GREY BROWN MUSHROOM



Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

 = Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

(US) 800-800-8586 | (CAN) 888-570-1113

pridemobility.com | prideliftchairs.com |    