

Specifications

Overall Weight	47 lbs (19Kg) (38 lbs with battery removed)
Dimension L x W x H	31 (785) x 25 (645) x 35 (889) - in. (mm)
Folded L x W x H	36 (921) x 25 (645) x 10 (254) - in. (mm)
Seat Bottom D x W x H	15 (380) x 16.5 (419) x 19 (485) measured from floor
Maximum User Weight	300 lbs (136Kg)
Speed Low Range	0-1 Mph (1.6 Km/h)
Speed Medium Range	0-3 Mph (5 Km/h)
Speed High Range	High 0-6 Mph (9.5 Km/h)
Reverse	0-1 Mph (1.6 Km/h)
Range	8 Miles (13 Km) flat continuous course, 5-6 miles typical outdoor use
Maximum Incline	10 Degrees (ADA compliant wheelchair ramps are 5° max)
Turning Radius	24 in. (610 mm)
Motor	Brushed DC 36V, 250W, 7A
Battery	36V, 6.6 AH, 240 watt-hour, Lithium-Ion 4 lbs(1.8Kg)
Charger	2.1A, Automatic, input 100-240V, charge time 4 approx. hours or less
Steering/Braking System	Two Lever Control System (tire brake and roller drive) requires both hands
Frame	6061 Aluminum and Nylon Composite
Front Wheels	6 inch non-marking PU solid
Rear Wheels	9 inch non-marking pneumatic (40PSI / 2.8Bar / 275kPa) easy change rims
Anti-Tip Wheels	38 inches (955mm) x 13 inches (330mm) x 28.5 inches (724mm) - 59 lbs (27 kgs)
Carton Size and Weight	38 inches (955mm) x 13 inches (330mm) x 28.5 inches (724mm) - 59 lbs (27 kgs)