

EVEN EASIER TO LOAD

We redesigned our best-selling model and broke it down into 3 easy to load pieces, all capable of fitting in the back of any vehicle. No need for specialty transport vehicles, vans, or expensive lifts.



CITY 2 PLUS Foldable Power Chair



Performance Specifications

Driving Range	15.5 mi / 25 km	Turning Radius	26"
Maximum Speed	4.5 mph / 7 kph	Stopping Distance	5 ft (at full speed)
Drive System	Two Motor Drive System	Obstacle Clearance	2"
Motor Type	250 Watt Brushless	Braking System	Electromagnetic
Incline Capability	12.0° / 21.3%	Max Weight Capacity	395 lbs

Dimensional Specifications

Unfolded Dimensions	44" L x 23" W x 40" H	Folded Dimensions	33" L x 23" W x 13" H
Chair Gross Weight	65 lbs	Front Wheel Diameter	8"
Chair Net Weight	58 lbs (no batteries)	Rear Wheel Diameter	12"
Seating Width	19.5" / 23"	Floor to Seat Height	19.5"
Seat Cushion Size	17" W x 18" D x 3" H	Back Support Height	20.5"

Other Specifications

Battery Type	Lithium-ion (24V 6Ah)	Charging Time	5 to 7 hours
Batteries Provided	Two (2)	Charger Input	100-240V, 50/60Hz ^(international)
Battery Watt Hour	144 Wh (Airline Approved)	Charger Output	29.4V 2.0A
Battery Charge Cycles	~1,000	Lap Belt	Provided
Tires	Non-Puncture	Horn	Built-In



Live Freely with Travel Buggy Electric Wheelchairs

Convenient to use. Quick to fold. Easy to transport.



1-855-56-BUGGY (28449)
us.travelbuggy.com
contact@travelbuggy.com



Your New Best Friend for Personal Mobility

You've waited long enough. It is time to see your friends and family, explore new places, maybe even travel to far off destinations, and it is all possible in our new CITY 2 PLUS. We've made our Best Selling CITY model even more convenient for you or your loved ones to fold, load, and go. It is time to get out there again. It is time to **join the Travel Buggy family.**

Easy to Fold and Transport

Folds quickly like a lawn chair to the size of a suitcase for convenience transport in any vehicle or for storage.



1-Step Fold Sequence

Net weight
58 lbs

Fits any vehicle

Will Travel!

Travel Buggy electric wheelchairs do not require any additional travel batteries or approvals like other "Travel Chairs". Our LiPo batteries are rated for safe airline travel by TCCA (CA), FAA (US) and various other aviation authorities.



Intelligently Designed

Our Travel Buggy electric wheelchairs are intuitive to drive and convenient to manage.



Multi-position reclining back rest (up to 135°) to relieve back pain & ensure comfortable all-day use



Batteries are easy to remove for airline travel



Entry and exit is made easier by a fold-away footrest & raiseable arms for optional side entry

Drive with Confidence Indoors and Outdoors

Do more indoors and outdoors with our CITY 2 PLUS model. Designed to handle what the outdoors can bring plus works great in tight spaces.



Use Indoors & Outdoors

4.5
mph

Top
Speed



15.5 MI
Long Driving
Range



Make Your Travel Buggy Unique To You

We offer more accessories than anyone else. From essentials like phone holders & extra storage to rehab accessories like head rests and leg extensions. We can accessorize your CITY 2 PLUS power chair to your desire.



The information contained herein is correct as at the time of printing, we reserve the right to alter specifications without prior notice