



FORCEMECH

FULL CARBON FIBER FRAME SERIES

lightest in the world!

Ultralite G10 Spec Sheet

Unfolded Dimensions (LWH)	36"x24"x35"
Folded Dimensions (LWH)	29"x12"x28"
Total Weight (w/o Battery)	26 Lbs
Battery Weight	2 Lbs
Turning Radius	28"
Weight Capacity	300 Lbs
Driving Distance	10 Miles
Max Climbing Angle	10°
Max Obstacle Height	1"
Max Speed	4 Mph
Minimum Breaking Distance (from max speed)	≤ 39"
Effective Seat Depth	16.5"
Seat Width	18"
Armrest Height	9.5"
Lithium Battery Specification	24V/10Ah
Ground Clearance	1.25"
Diameter of Front Wheel	6.5"
Diameter of Rear Wheel	11"
Motor Specification	24V 200W



Scan for
Complete Info!

Mobility FREEDOM at Your Finger Tips!