



# FORCEMECH

FULL CARBON FIBER FRAME SERIES

**Vacation Ready!**

## Carbon F1 Spec Sheet

Unfolded Dimensions (LWH)	27"x24"x38"
Folded Dimensions (LWH)	28"x13"x24"
Total Weight (w/o Battery)	33 Lbs
Battery Weight	2 Lbs
Turning Radius	24"
Weight Capacity	300 Lbs
Driving Distance	10 Miles
Max Climbing Angle	9°
Max Obstacle Height	1"
Max Speed	4 Mph
Minimum Breaking Distance (from max speed)	≤ 39"
Effective Seat Depth	17"
Effective Seat Width	20"
Armrest Height	9"
Lithium Battery Specification	24V/10Ah
Ground Clearance	1.25"
Diameter of Front Wheel	7"
Diameter of Rear Wheel	8.5"
Motor Specification	24V 200W



1-877-903-6723

WWW.FORCEMECH.COM

Forcemech Carbon F1 has the Smallest Folded  
Signature of any Powerchair in the World!