

**icare**  
Medical Group

**PRODUCT  
CATALOG**  
NORTH AMERICA



# LIVING WELL

*Transforming the daily lives of people in need.*



## Wellbeing Champions

Of our people, end users,  
carers, healthcare  
professionals and dealers.



## Purposeful Design

End user appeal with  
clinical performance.



## Collaboration Care

We speak to people to  
find solutions that work.



## Authentic Heart

Real passion and care,  
humble and honest team,  
part of your family.



# GLOBAL REACH




## Distribution

With dealers and service agents across North America, New Zealand and Australia we are delivering solutions to homes & facilities everyday. Our distribution centers in Nebraska and California allow us to quickly supply to anywhere in North America.

**ICARE MEDICAL GROUP**  
[icaremedicalgroup.com](http://icaremedicalgroup.com)  
[us.support@icaremedicalgroup.com](mailto:us.support@icaremedicalgroup.com)


+1 308 351 3000

## PRODUCTS



5

BEDS & PARTNER  
BEDS




15

BED  
ACCESSORIES




23

MATTRESSES



30

PILLOWS &  
POSITIONING



34

MATTRESS  
PROTECTION  
& SHEETING



39

RESOURCES



# BEDS & PARTNER BEDS

---




# ICARE IC111

The IC111 Bed is our low risk extra mobility option. The IC111 is a static height bed with head lift and knee break. This allows the user to sit up and provides support through the legs.

Available in a range of sizes and upholstery options to suit any bedroom or room setting.



## BED CODES & SIZES




### TWIN XL

Code: IC111TTL

Mattress Surface Size: 80" x 35"

Base Size: 85" x 38"



### DUAL KING

For a dual king order 2x IC111TTL

Mattress Surface Size: 80" x 35"

Base Size: 85" x 38"

- Individual hand control per side
- Requires two Twin XL headboards

## KEY FUNCTIONS



Head Lift



Knee Break

## COLOR CHOICES




ONYX



STONE



SMOKE



SILVER



VINYL EBONY  
(wipeable)



VINYL DOVE  
(wipeable)



## KEY FEATURES



440LBS Safe  
Working Load per  
side



Fixed Base Height  
9 Inches



Available with  
Headboards and  
Footboards



External Rechargeable  
Battery Back-up



Lockable Heavy Duty  
Casters



7 Inches Underbed  
Clearance



End and Corner  
Mattress Retainer Bars



5 Year Warranty



Backlit Wipeable  
Hand Piece



No gap between  
mattresses on Queen  
& King Beds

## TECHNICAL SPECIFICATIONS

- Safe working load - The safe working load is 400LBS per side and includes any model of Icare mattress.
- Hand Piece - Backlit for night time use. Easy press buttons. Wired stretch cord. Hook attachment for hanging to side of bed. Wipeable.
- Fabric - Hard wearing commercial grade fabric. Vinyl is treated with antibacterial properties.
- Mattress Surface - Strong metal mesh under mattress allowing increased airflow.
- Actuators - 1x 6000 newton - 24 watt head adjust actuator. 1x 6000 newton - 24 watt knee adjust actuator.
- No pull motors for added safety.
- Wheels - large lockable heavy duty swivel casters on each corner.
- Head Lift angle 62°
- Knee break angle 42°
- Fixed base height - 9" from floor to mattress platform.

## COMPATIBLE ACCESSORIES AVAILABLE

- ✓ High Side Rail
- ✓ Low Side Rail
- ✓ U-Assist Rail
- ✓ Headboard
- ✓ Footboard
- ✓ Slide-Out Side Rail



# ICARE IC222

The IC222 bed combines style with the essential features of a homecare bed. This three function design is supported with heavy duty locking casters giving the bed 8" of under-bed clearance for lifting hoists.



## KEY FUNCTIONS



Head Lift




Knee Break




Height Adjustable  
12" - 24"

## BED CODES & SIZES



### TWIN XL

Code: IC222TXL  
Mattress Surface Size: 80" x 35"  
Base Size: 85" x 38"



### DUAL KING

For a dual king order 2x IC222TXL  
Mattress Surface Size: 80" x 35"  
Base Size: 85" x 38"  
• Individual hand control per side  
• Requires two Twin XL headboards

## COLOR CHOICES




ONYX



STONE



SMOKE



SILVER



VINYL EBONY  
(wipeable)



VINYL DOVE  
(wipeable)



## KEY FEATURES



390LBS Safe Working Load



Strong Breathable Mattress Platform



Available with Headboards and Footboards



External Rechargeable Battery Back-up



Lockable Heavy Duty Castors



End and Corner Mattress Retainer Bars



5 Year Warranty



Backlit Wipeable Hand Piece

## TECHNICAL SPECIFICATIONS

- Safe working load - 390LBS. The safe working load includes any model of Icare mattress.
- Hand Piece - Backlit for night time use. Easy press buttons. Wired stretch cord. Hook attachment for hanging to side of bed.
- Wipeable.
- Fabric - Hard wearing commercial grade fabric. Vinyl is treated with antibacterial properties.
- Mattress Surface - Strong metal mesh under mattress allowing increased airflow
- No pull motors for added safety.
- Wheels - 4x heavy duty lockable swivel casters.
- Head lift angle 62° and knee break angle 42°
- Max height - 24" from floor to top of mattress platform.
- Min height - 12" from floor to top of mattress platform.

## COMPATIBLE ACCESSORIES AVAILABLE


- ✓ High Side Rail
- ✓ Low Side Rail
- ✓ U-Assist Rail
- ✓ Headboard
- ✓ Footboard
- ✓ Slide-Out Side Rail

# ICARE IC333

The IC333 is a four function bed in a range of sizes, colors and accessories. This innovative design combines aesthetic appeal with all the essential features. With its ultra-low height, four-functions, and ample under-bed clearance, this bed offers the perfect blend of aesthetics and functionality .



## BED CODES & SIZES



### TWIN XL

Code: IC333TXL  
Mattress Surface Size: 80" x 35"  
Base Size: 85" x 38"




### FULL XL

Code: IC333FXL  
Mattress Surface Size: 80" x 53"  
Base Size: 85" x 55"  
• One hand control



### QUEEN

Code: IC333Q  
Mattress Surface Size: 80" x 60"  
Base Size: 85" x 62"  
• One hand control



### DUAL KING

For a dual king order 2x IC333TXL  
Mattress Surface Size: 80" x 35"  
Base Size: 85" x 76" (as partner)

## KEY FUNCTIONS



Head Lift



Knee Break



Height Adjustable  
9" - 26"



Trendelenburg and  
Reverse Trendelenburg

## COLOR CHOICES



ONYX



STONE



SMOKE



SILVER



VINYL EBONY  
(wipeable)



VINYL DOVE  
(wipeable)



## KEY FEATURES



660LBS SWL  
Full XL and Queen



440LBS SWL  
Twin XL



Unlimited Under Bed  
Clearance



Available with  
Headboards and  
Footboards



External Rechargeable  
Battery Back-up



Lockable Heavy Duty  
Casters



Ultra Low Design  
Reducing Risk of Falls



End and Corner  
Mattress Retainer Bars



10 Year Warranty



Backlit Wipeable  
Hand Piece



No gap between  
mattresses Dual King  
Beds



Strong Breathable  
Mattress Platform



Optional Three Sided  
Surround for  
Wheelchair Transfers

## TECHNICAL SPECIFICATIONS

- Safe working load - 440LBS on long single and king single sizes and 660LBS for double and queen sizes.
- The safe working load includes any model of Icare mattress.
- Hand Piece - Backlit for night time use. Easy press buttons. Wired stretch cord. Hook attachment for hanging to side of bed. Wipeable.
- Fabric - Hard wearing commercial grade fabric. Vinyl is treated with antibacterial properties.
- Mattress Surface - Strong metal mesh under mattress allowing increased airflow
- Actuators - 1x 6000 newton - 24 watt head adjust actuator. 1x 6000 newton - 24 watt knee adjust actuator. 2x 6000 newton - 24 watt height adjust actuators (2x 8000 newton on Double and Queen sizes)
- No pull motors for added safety.
- Wheels - 8x heavy duty swivel castors with one lockable castor on each corner.  
Head lift angle 62° and knee break angle 42°
- Max height - 26" from floor to top of mattress platform.
- Min height - 9" from floor to top of mattress platform.

NOTE: Always remember that TWO Twin XL beds together make a Dual King size bed.

## COMPATIBLE ACCESSORIES AVAILABLE

- ✓ High Side Rail
- ✓ Low Side Rail
- ✓ U-Assist Rail
- ✓ Headboard
- ✓ Footboard
- ✓ Slide-Out Side Rail



# ICARE IC100 COMPANION BASE

The IC100 Bed is a heavy duty bed, with a static mattress platform. This bed provides a strong base for a partner or parent to sleep alongside the client. Designed for clinical applications with medical grade casters and a wipeable mattress platform.



440LBS Weight Rating



Lockable Heavy Duty Casters

## TECHNICAL SPECIFICATIONS

- Weight rating – The max weight limit is 440LBS, includes any model of Icare mattresses.
- Fabric – Hard wearing commercial grade fabric. Vinyl is treated with antibacterial properties.
- Mattress Surface – Strong metal mesh under mattress allowing increased airflow.
- Wheels – large lockable heavy duty swivel casters on each corner.
- Fixed base height – 14" from floor to mattress platform.
- Under Bed clearance – 7"




Strong Breathable Mattress Platform



10 Year Warranty

## BED CODES & SIZES



### TWIN XL

Code: IC100TXL

Mattress Surface Size: 80" x 35"

Base Size: 85" x 38"

## COLOR CHOICES



ONYX



STONE



SMOKE



SILVER



VINYL EBONY  
(wipeable)



VINYL DOVE  
(wipeable)



# ICARE PARTNER SET-UPS


All Icare Beds can be made for couples to sleep side by side with no gap between mattresses. You simply choose in 3 steps.

1. Choose which model each partner wants.
2. Choose the bed size
3. Choose which side each partner will be sleeping on. This can be filled out on the Icare Partner Bed Assessment Script Form.


Note: Headboards are sold as a separate item and can only be sold to fit each side.



## SIDE X SIDE SIZE CONFIGURATIONS:



Note: This makes a Dual king size bed



This form can be completed using our online form, or scanned in using PDF document.



# ICARE PARTNER SET-UPS



Center rails are upholstered and recessed.



Side x Side configurations.



No gap between the mattresses



Separate headboards for each bed.




Mix and match bed models.



Left or right side sleeping preferences.

## MODEL CONFIGURATION OPTIONS:



- No gap between the mattresses
- Fully upholstered surround
- Separate headboards
- Compatible with all side rails and accessories





# BED ACCESSORIES

---



## U-ASSIST SIDE RAIL

A unique design rail, slim line designed for bed users to grab anywhere. Requested by Occupational Therapists as a variation to the Bed Stick. Secure mount to the bed, multiple grip points and assists bed users with transfers.

Charcoal color, slip on attach system. Compatible with any Icare Beds. Position anywhere along side of bed as needed.

Width 6"

Height above mattress platform 21"

Max user weight 440LBS (static load maximum 56 pounds)



Compatible with:

- IC100
- IC111
- IC222
- IC333

PRODUCT CODE	
Description Code	Icare U-Assist Side Rail ACUSR

## EXTENDABLE SIDE RAIL

The Icare Extendable Side Rail can be mounted on either side of the bed.

Charcoal colour, slip on attach system with four rails for maximum safety. Position anywhere along side of bed as needed. Locks in two positions, standard width or fully extended.

Width 27" - 37"

Height above mattress platform 17"

Space between rails 4"

Max user weight 440LBS



Compatible with:

- IC100
- IC111
- IC222
- IC333

PRODUCT CODE	
Description Code	Icare Extendable Side Rail ACESR



## LOW SIDE RAIL

Provides the user support and stability, preventing bed falls.

Charcoal color, slip on attach system with three rails for maximum safety. Position anywhere along side of bed as needed.

Width 27"

Height above mattress platform 17"

Space between rails 4.5"

Max user weight 440LBS (static load maximum 56 pounds)



Compatible with:

- IC100
- IC111
- IC222
- IC333

PRODUCT CODE	
Description Code	Icare Low Side Rail ACLSR

## HIGH SIDE RAIL

Provides the user support and stability, especially for getting in and out of bed.

Charcoal color, slip on attach system with four rails for maximum safety. Compatible with any Icare Beds. Position anywhere along side of bed as needed.

Width 13.75"

Height above mattress platform 23"

Space between rails 4.5"

Max user weight 440LBS (static load maximum 56 pounds)



Compatible with:

- IC100
- IC111
- IC222
- IC333

PRODUCT CODE	
Description Code	Icare High Side Rail ACHSR

## FULL LENGTH SIDE RAIL

Full length fold down rails available for the left and right hand side. Release is designed for caretakers use.

Charcoal color, slip on attach system with four rails for maximum safety. Left and right sides available. Rail folds towards the foot end of the bed.

Width 54.5"

Height above mattress platform 19.5"

Space between rails 3.25"

Max user weight 440LBS (static load maximum 56LBF)

### PRODUCT CODES

Description Code	Icare Full Length Rail - Right ACFSRR
Description Code	Icare Full Length Rail - Left ACFSRL
Description Code	Icare Full Length Rail - Pair ACFSRP



Prescribe rail from clients perspective:

Rail to be on clients left = left rail

Rail to be on clients right = right rail.

Compatible with:

- IC100
- IC111
- IC222
- IC333

## OVER BED POLE

Over Bed Pole to help with turning, patient transfer and bed care. Free standing, stable design with heavy duty steel and strong bracing brackets. Adjustable handle height and anti-slip strap cradle.

Charcoal color. Compatible with any Icare Beds. Low profile base to fit neatly under ultra low beds, suitable for all bed sizes. Minimal assembly required.

Height from floor to top of rail 74.75"

Floor Foot Print: H design

- Width under bed 29.5"
- Length under bed 35.5"

Rail at top from back of stand to tip of rail – 28"

Maximum weight applied 220LBS

### PRODUCT CODE

Description Code	Icare Pull Up Bar ACPUB
---------------------	----------------------------



Compatible with:

- IC100
- IC111
- IC222
- IC333



# PADDED SIDE RAIL COVERS

Icare padded rail covers for protection on the Icare side rails. The Icare padded covers are covered with Zerotec wipeable medical fabric.

Zerotec super-soft two way stretch reducing sheer and skin damage. This medical fabric combines patient care and technology to provide unparalleled comfort.

Completely wipeable, full rail height padded sleeve secured in place using velcro.



## PRODUCT CODES

Description Code	Icare Bed Stick/ U-Assist Cover ACSRCB Fits rail size: 4" x 20.25"
Description Code	Icare Low Side Cover ACSRCL Fits rail size: 27" x 16.5"
Description Code	Icare High Side Cover ACSRCH Fits rail size: 13.75" x 23"
Description Code	Icare Full Length Cover ACSRCFL Fits rail size: 54.5" x 19.25"

# BLANKET CRADLE

Blanket cradles are designed to support the weight of bed linen. This can relieve the lower legs and feet from pressure while also keeping sheets or blankets away from sores or wounds. Blanket cradles can also promote air circulation to help keep skin dry. Easy height adjustment system, and removable. Charcoal color, easy to mount to the bed base/mattress platform, moves with the mattress platform.

Measured from mattress platform 17" to 19.75".  
Width 16.5"



## PRODUCT CODE

Description Code	Icare Blanket Cradle ACOBC
---------------------	-------------------------------

Compatible with:

- IC100
- IC111
- IC222
- IC333

## HEADBOARDS

Headboards complete the look in the bedroom with matching colors and fabrics. Headboards can prevent the pillow or other bedding from falling off the head of the bed when the bed is functioning.

Mounted to the bed using 4 hand turn bolts.  
Total height: 27.5" (20" above mattress platform)

Available in all fabric colors.

### PRODUCT CODES

Description Code	Icare Headboard - Twin XL HB333TXL
Description Code	Icare Headboard - Full XL HB333FXL
Description Code	Icare Headboard - Queen HB333Q



Compatible with:

- Companion
- IC100
- IC111
- IC222
- IC333

## FOOTBOARDS

Footboards can prevent the mattress and bed linen from falling off the bed or slipping out of reach. They also provide a soft barrier when carers or users move around the bed room. Easy to mount using metal/plastic bracket system. Not designed to bear weight, do not use as a side rail.

Total height 20" (12" above mattress platform)

Available in all fabric colors.

### PRODUCT CODES

Description Code	Icare Footboard - Twin XL FB333TXL	Description Code	Icare Footboard - Full XL FB333FXL
Description Code	Icare Footboard - Queen FB333Q		

Compatible with:

- Companion
- IC100
- IC111
- IC222
- IC333





# OVER BED TABLE

The Icare Over Bed Table provides a solid foundation for work, recreation or meals, conveniently sits over any bed and mattress. Can also be used over a chair. Easy care wipe down split table top, providing large tilting workspace and separate 'side table' all in one. The tilting table top is reversible.

Height Adjustment Range	28.75" to 43"
Tilt Top Surface Area	18.5" x 12"
Tilt Top Stepped Range	0 - 70°
Fixed Side Table Measurements	18.5" x 6"
Table Top Load Weight Limit	65 LBS


Standard Height (larger wheels)	3.5"
Low-Profile Height	1.75"



Reversible tilting table top.



Easy to assemble.



Gas assisted height adjustable.



Universal Compatibility

## PRODUCT CODES

Description Code	Icare Over Bed Table - Standard ACOBTH
------------------	--



## PRODUCT CODES

Description Code	Icare Over Bed Table - Low Profile ACOBTL
------------------	---



## OVER BED TABLE C SHAPE BASE

The Icare Over Bed Table C Shape Base, fits easier over a chair and can be used with other homedcare beds. Easy care wipe down split table top, providing large tilting workspace and separate 'side table' all in one. The tilting table top is reversible.

C-Shape Base Height 3.5"  
Base Width 31.5"



### PRODUCT CODES

Description Code	Icare Over Bed Table - C Shape ACOBTC
------------------	--



## ORGANIZER TRAY

Soft durable silicone, manufactured to provide a non-slip surface, slips over the edge rails of the Icare Over Bed Table. Compatible with Icare Over Bed Tables.

Four compartments, designed to securely hold personal items in place.



### PRODUCT CODE

Description Code	Icare Organizer Tray ACOT
------------------	------------------------------





# ICARE PRESSURE CARE MATTRESSES

---



# IC15 ACTIVEX™ MATTRESS


Suitable for the prevention and treatment of pressure ulcers, the IC15 offers firm comfort and support for ease of mobility.

RISK CATEGORY:

**HIGH RISK**

COMFORT CATEGORY:

**FIRM FEEL**



## KEY FEATURES



PRESSURE  
REDUCING



TEMPERATURE  
REGULATING




BREATHABLE



ANTI SHEAR



MOISTURE  
WICKING




WITHOUT GOING TOO FAR

XP Support prevents you from bottoming out the mattress even on adjustable beds.

IMMERSE YOURSELF....

ActiveX® provides optimal immersion and envelopment for long term comfort and pressure relief.



DONT SWEAT IT

Bluetec® and Vant® fabric technologies combine to keep your skin dry and cool.




## MORE FEATURES



10 Year  
Warranty



660LBS  
Weight Rating



6" Mattress  
Height



Compatible with  
Zerotec Waterproof  
Covers



Flip Free Design



Zip-off Washable  
Covers



Multi Layer  
Immersion



Anti Slip Base  
Layer




Designed for use on  
adjustable beds




20+ Waterlow  
Rating

## BED CODES &amp; SIZES



**TWIN XL**  
Code: IC15TXL  
Mattress Size: 80" x 35"



**FULL XL**  
Code: IC15FXL  
Mattress Size: 80" x 35"



**QUEEN**  
Code: IC15Q  
Mattress Size: 80" x 60"

## SPECIFICATIONS

ActiveX™ pressure relieving material  
Vant ventilated sides and base fabric  
Bluetec temperature regulating cover  
Option Zerotec waterproof cover (see pg. 33)  
Flip free design  
Zip off machine washable covers  
100 - 660LBS user weight  
Standard cover MVTR: NA.  
Medical Cover MVTR (g/24h·m2): 1635.0 Average  
Medical Cover BS7175: Crib 5  
Waterlow ≤20

Feel Rating	Firm
Mattress Height	150mm
Risk Category	High Risk
User Weight	100LBS - 660LBS
Warranty	10 Year Warranty

## PRESSURE CARE

Risk Category High

This rating is based on our average peak pressure of <24mmHg. For some high risk clients Icare Zerotec covers are recommended.

Icare risk category ratings are based on average peak pressure, Therapist recommendation, case studies, industry comparison and the Waterlow scale. For some high risk clients Icare Zerotec covers are recommended. The best way to assess a client on a pressure reducing surface is to monitor pressure points with visual checks over a period of 5 days.

## CARE &amp; CLEANING

To launder the cover in a washing machine, set the main wash to a minimum of 150°F for no less than 10 minutes. The cover can be tumble dried on low heat and the temperature must not exceed 140°F. The foam can be autoclaved at 273°F.

# IC20 ACTIVEX™ MATTRESS


Balancing pressure care and comfort, the IC20 is a medium feel, providing excellent immersion, support and in-bed mobility.

RISK CATEGORY:

HIGH RISK

COMFORT CATEGORY:

MEDIUM FEEL



## KEY FEATURES



PRESSURE  
REDUCING



TEMPERATURE  
REGULATING




BREATHABLE



ANTI SHEAR



MOISTURE  
WICKING




## SPRING FREE DESIGN

Adjustable bed compatibility.  
Low average peak pressure.  
No push Back.

## BODY HEAT PRESSURE CARE

Activex™ uses body heat to soften under hot areas, preventing pressure ulcers.



## COOLING TECHNOLOGY

Bluetec® and Vant® fabric technologies combine to keep your skin dry and cool.




## MORE FEATURES



10 Year  
Warranty



660LBS  
Weight Rating



8" Mattress  
Height



Compatible with  
Zerotec Waterproof  
Covers



Flip Free Design



Zip-off Washable  
Covers



Multi Layer  
Immersion



Anti Slip Base  
Layer




Designed for use on  
adjustable beds




20+ Waterlow  
Rating


## BED CODES &amp; SIZES



**TWIN XL**  
Code: IC2TXL  
Mattress Size: 80" x 35"



**FULL XL**  
Code: IC2FXL  
Mattress Size: 80" x 53"



**QUEEN**  
Code: IC2Q  
Mattress Size: 80" x 60"

## SPECIFICATIONS

ActiveX™ pressure relieving material  
Vant ventilated sides and base fabric  
Bluetec temperature regulating cover  
Option Zerotec waterproof cover (see pg. 33)  
Flip free design  
Zip off machine washable covers  
66LBS - 440LBS user weight  
Standard cover MVTR: NA.  
Medical Cover MVTR (g/24h-m<sup>2</sup>): 1635.0 Average  
Medical Cover BS7175: Crib 5  
Waterlow ≤20

Feel Rating	Medium
Mattress Height	200mm
Risk Category	High Risk
User Weight	66LBS - 440LBS
Warranty	10 Year Warranty

## PRESSURE CARE

Risk Category High

This rating is based on our average peak pressure of <19mmHg. For some high risk clients Icare Zerotec covers are recommended.

Icare risk category ratings are based on average peak pressure, Therapist recommendation, case studies, industry comparison and the Waterlow scale. For some high risk clients Icare Zerotec covers are recommended. The best way to assess a client on a pressure reducing surface is to monitor pressure points with visual checks over a period of 5 days.

## CARE &amp; CLEANING

To launder the cover in a washing machine, set the main wash to a minimum of 150°F for no less than 10 minutes. The cover can be tumble dried on low heat and the temperature must not exceed 140°F. The foam can be autoclaved at 273°F.

# IC30 ACTIVEX™ MATTRESS

The IC30 is our all in one mattress. Incorporating pressure care, bed mobility and most importantly a very comfortable sleeping surface. Highly effective in the treatment and prevention of pressure sores.

RISK CATEGORY:

**HIGH RISK**

COMFORT CATEGORY:

**MEDIUM FEEL**



Aged Care




Every Day Care



Disability




Facilities



Cooling  
Comfort  
Care

## TRANSFERS

Transferring to and from the IC30 mattress is assisted with the XP Support mattress edge and anti-hear BlueTec mattress surface cover.



## SLEEP QUALITY

Quality sleep is an important factor to wellbeing and healing. The IC30 mattress creates an ideal balance for medical and pressure care needs, sleep quality and longevity.



## PRESSURE CARE

The IC30 Mattress is designed to prevent and reduce the risk of pressure injuries. The ActiveX™ sleep surface provides partial immersion with envelopment, creating an excellent average peak pressure, while maintaining bed mobility and lateral movement. The mattress microclimate is regulated by using BlueTec, gel infusion and charcoal.



## FEATURES



Pressure Relieving

ZeroTec Covers  
availableActiveX™ Gel  
temperature regulating

Anti-shear Materials

Moisture Wicking &  
Breathable FabricActiveX™ Charcoal  
reduces odors and  
is microbial.Multi Layer  
ImmersionReduces Rollout  
RiskContour-Latch™  
attaches to most  
profiling beds.20+ Waterlow  
Rating

## BED CODES &amp; SIZES

TWIN XL  
Code: IC30TXL  
Mattress Size: 80" x 35"FULL XL  
Code: IC30FXL  
Mattress Size: 80" x 53"QUEEN  
Code: IC30Q  
Mattress Size: 80" x 60"

## SPECIFICATIONS

ActiveX™ pressure relieving material  
 Vant ventilated sides and base fabric  
 Bluetec temperature regulating cover  
 Option Zerotec waterproof cover (see pg. 33)  
 Flip free design  
 Zip off machine washable covers  
 66LBS - 440LBS user weight  
 Standard cover MVTR: NA.  
 Medical Cover MVTR (g/24h-m<sup>2</sup>): 1635.0 Average  
 Medical Cover BS7175: Crib 5  
 Waterlow ≤20

Feel Rating	Medium
Mattress Height	200mm
Risk Category	High Risk
User Weight	45LBS - 502LBS
Warranty	10 Year Warranty

## PRESSURE CARE

Risk Category High

This rating is based on our average peak pressure of <19mmHg. For some high risk clients Icare Zerotec covers are recommended.

Icare risk category ratings are based on average peak pressure, Therapist recommendation, case studies, industry comparison and the Waterlow scale. For some high risk clients Icare Zerotec covers are recommended. The best way to assess a client on a pressure reducing surface is to monitor pressure points with visual checks over a period of 5 days.

## CARE &amp; CLEANING

To launder the cover in a washing machine, set the main wash to a minimum of 150°F for no less than 10 minutes. The cover can be tumble dried on low heat and the temperature must not exceed 140°F. The foam can be autoclaved at 273°C.

# PILLOWS & POSITIONING

---





# CONFORM PILLOW

The Icare Conform Pillow is an Adjustable Pillow designed to provide multiple options with maximum comfort and ultimate positioning.



## MULTIPLE LAYERS & HEIGHTS

The Conform Pillow features the ActiveX material which makes up the modules of this pillow.



## HYPERSOFT COVER

The cover features soft and gentle textures, hypoallergenic materials, and is designed to reduce noise and tactile stimulation.



## CHARCOAL INFUSED ACTIVEX™

Charcoal material is durable and known to be antimicrobial. It can assist with heat dispersion, drawing away moisture, and odor resistant.



## UNIQUE FIBER FILLED ENCASEMENT

The unique hollow core fibre encased cover gives the user a soft breathable comfort layer.



## PRESSURE REDUCING

ActiveX is a high grade type of elastic foam which is responsive to temperature around inflamed areas.

## ONE PILLOW, MANY SOLUTIONS.

Made up of multiple layers this pillow is truly adjustable (the first of its kind). The Conform provides complete customisation to suit all sleepers and sleeping positions, and assists with positioning, alignment, and support as required.



### PRODUCT CODE

Description Code	Icare Conform Pillow ICP6
------------------	---------------------------

## CLASSIC PILLOW

The Icare Classic pillow is made up of a high density ActiveX core in a traditional pillow shape, a cooling gel insert and a soft puff quilt style cotton cover. Zip off washable cover. Suitable for back and side sleepers.

The Icare Classic Pillow is made with an Elastex material cover. This fabric has a high elasticity in every direction enabling it to conform with head and neck contour.




Size: 23.6" (L) x 16.5" (W) x 4.3" (H)

PRODUCT CODE	
Description Code	Icare Classic Pillow ICP1



### PRESSURE REDISTRIBUTION

The Classic Pillow features the ActiveX material which makes up the core of this pillow. ActiveX is a high grade type of elastic foam which is responsive to temperature around inflamed areas.




### COOLING GEL INSERT

Under the fabric this pillow features a cooling gel insert. This helps to regulate temperature and will conform with head and neck contour.

## CONTOUR PILLOW

The Icare Contour neck pillow is made up of a high density ActiveX core in a contoured pillow shape, a cooling gel insert and a soft puff quilt style cotton cover. This pillow is recommended for back and neck pain sufferers. Higher side for back and side sleepers and lower side for front sleepers.

The Icare Contour Pillow is made with an Elastex material cover. This fabric has a high elasticity in every direction enabling it to conform with head and neck contour and is zipped to allow for washing.




Size: 23.6" (L) x 16.5" (W) x 3.5"-4.3" (H)

PRODUCT CODE	
Description Code	Icare Contour Pillow ICP2



### PRESSURE REDISTRIBUTION

The Contour Pillow features the ActiveX material which makes up the core of this pillow. ActiveX is a high grade type of elastic foam which is responsive to temperature around inflamed areas.



### COOLING GEL INSERT

Under the fabric this pillow features a cooling gel insert. This helps to regulate temperature and will conform with head and neck contour.



## CLOUD PILLOW

The Icare Cloud pillow is made up of a high density ActiveX ventilated core with a microfibre filled cotton cover. The ActiveX core is removable to make the Cloud pillow even softer. Zip off washable cover. Suitable for all sleepers.

The Cloud pillow features a unique microfibre filled cover. This gives an extremely soft feel with a high support and pressure relieving core. The core and the cover feature pinhole ventilation and hollow fibres.



Size: 24.8" (L) x 16.5" (W) x 5" (H)

### PRODUCT CODE

Description Code	Icare Cloud Pillow ICP3
------------------	-------------------------



### PRESSURE REDISTRIBUTION

The Cloud Pillow features ventilated ActiveX material which makes up the core of this pillow. ActiveX is a high grade type of elastic foam which is responsive to temperature around inflamed areas.



### REMOVABLE INSERT

The pinhole ventilated ActiveX core is covered with an inner fabric which is able to be removed to make the Cloud Pillow even softer.

## CURVE PILLOW

The Icare Curve neck pillow is made up of a high density ActiveX core in an ergonomic shape designed to support the shoulders and neck. This pillow is recommended for back and neck pain sufferers. Very low profile, suitable for children, petite clients and those on bed rest.

The Icare Curve Pillow is made with an Elastex material cover. This fabric has a high elasticity in every direction enabling it to conform with head and neck contour.



Size: 22" (L) x 13.4" (W) x 2"-3" (H)

### PRODUCT CODE

Description Code	Icare Curve Pillow ICP4
------------------	-------------------------



### PRESSURE REDISTRIBUTION

The Curve Pillow features ventilated ActiveX material which makes up the core of this pillow. ActiveX is a high grade type of elastic foam which is responsive to temperature around inflamed areas.



### VENTILATED SIDE FABRIC

Vant fabric is a three-dimensional polypropylene material that has high air permeability and is antibacterial. Using Vant on the sides of the pillow greatly increases air movement.

# MATTRESS PROTECTION AND SHEETING

---





# MATTRESS PROTECTOR

The Icare waterproof Mattress Protector has a Tencel surface with a noiseless breathable PU backing. The microporous backing provides a breathable yet waterproof layer allowing moisture to dissipate. Designed in a fitted sheet style the Icare protectors provide protection for the top of the mattress from liquid and stains. 12.5" drop.



100% Waterproof top layer sheet



Home or facility use



Hygiene cover for general bed protection

## CARE & CLEANING

This item is machine washable on delicate cycle, line dry in the shade. We do not recommend machine drying, ironing or dry-cleaning. Ensure the item has fully dried before reusing.

## FITTED SHEET STYLE



## PRODUCT CODE

Description Code	Mattress Protector - Twin XL ICMPTXL 79" x 35"
Description Code	Mattress Protector - Full XL ICMPFXL 79" x 54"
Description Code	Mattress Protector - Queen ICMPQ 80" x 60"

## WITHOUT PROTECTION



## WITH MATTRESS PROTECTION



# MATTRESS COVER

The Icare fully enclosed cover has a Tencel surface with a noiseless breathable PU backing. The microporous backing provides a breathable yet waterproof layer allowing moisture to dissipate. The Mattress Cover fully encloses the mattress and zips up protecting the user and mattress from liquid and stains. 13.75" drop.

FULLY ENCLOSED WITH ZIP



100% Waterproof



Home or facility use



Hygiene cover for general bed protection



Prevents bed bugs & dustmites

## CARE & CLEANING

This item is machine washable on delicate cycle, line dry in the shade. We do not recommend machine drying, ironing or dry-cleaning. Ensure the item has fully dried before reusing.

### PRODUCT CODE

Description    Mattress Cover - Twin XL  
Code            ICMCTXL    79" x 35"

Description    Mattress Cover - Full XL  
Code            ICMCFXL    79" x 54"

Description    Mattress Cover - Queen  
Code            ICMCQ      80" x 60"

### WITHOUT PROTECTION



### WITH MATTRESS PROTECTION





# MEDICAL COVERS

The medical mattress cover is made from Zerotec fabric. A medical barrier fabric that has high elasticity in every direction which reduces shear on pressure redistribution surfaces. It is waterproof and wipeable but has microporous breathable properties. This cover has welded seams and waterfall flaps for fluid control.



Strong zipper



Waterfall flaps



## ZEROTEC COVER

Zerotec fabric is a medical barrier fabric that has high elasticity in every direction which reduces shear on pressure redistribution surfaces. It is waterproof and wipeable but has microporous breathable properties. This cover has welded seams and waterfall flaps for fluid control.



## SPECIFICATIONS

- Medical Cover MVTR (g/24h-m2): 1635.0 Average
- Medical Cover BS7175: Crib 5

## CARE & CLEANING

To launder the cover in a washing machine, set the main wash to a minimum of 150°F for no less than 10 minutes. The cover can be tumble dried on low heat and the temperature must not exceed 140°F.

### PRODUCT CODE

Description 6" Mattress Cover - Twin XL  
Code ICMC15TXL  
35" (W) x 80" (L) x 6" (D)

Description 6" Mattress Cover - Full XL  
Code ICMC15FXL  
53" (W) x 80" (L) x 6" (D)

Description 6" Mattress Cover - Queen  
Code ICMC15Q  
60" (W) x 80" (L) x 6" (D)

### PRODUCT CODE

Description 8" Mattress Cover - Twin XL  
Code ICMC20TXL  
35" (W) x 80" (L) x 8" (D)

Description 8" Mattress Cover - Full XL  
Code ICMC20FXL  
53" (W) x 80" (L) x 8" (D)

Description 8" Mattress Cover - Queen  
Code ICMC20Q  
60" (W) x 80" (L) x 8" (D)

## ADJUSTABLE BED SHEET SETS

Bed sheets that are made specially for Icare profiling beds in all sizes. Sheet Sets feature a 1000 thread count sateen which minimizes shear. Suitable for all adjustable beds, large 14" drop for a secure fit to prevent becoming untucked during use.



### CARE & CLEANING

Icare Sheets are machine washable, we recommend line dry in the shade. Suitable for machine drying, ironing or dry-cleaning. Ensure the item has fully dried before storing or reusing.

PRODUCT CODE				
Description Code	Twin XL Sheet Set	White ICSSTXL	Charcoal ICSSCTXL	Fits bed size: 35" (W) x 80" (L) 1x Flat Sheet, 1x Fitted Sheet, 1x Pillow Case
Description Code	Full XL Sheet Set	White ICSSFXL	Charcoal ICSSCFXL	Fits bed size: 53" (W) x 80" (L) 1x Flat Sheet, 1x Fitted Sheet, 2x Pillow Cases
Description Code	Queen Sheet Set	White ICSSQ	Charcoal ICSSCQ	Fits bed size: 60" (W) x 80" (L) 1x Flat Sheet, 1x Fitted Sheet, 2x Pillow Cases
Description Code	Dual King Sheet Set	White ICSSDK	Charcoal ICSSCDK	Fits bed size: 35" (W) x 80" (L) per side 1x Flat Sheet, 2x Fitted Sheets, 2x Pillow Cases



# RESOURCES





## DEALER INFORMATION

Icare Medical Group has a range of dealer resources available from images and point of sale, we can provide; Mattress sample cubes, brochures, in-store banners and posters, videos and images, user manuals and warranty cards, wall graphics, display mats and covers, product and fabric samples, stands and shelving.

To reorder marketing material please ask for an order form.



To become a dealer or request a meeting with a product specialist, visit our website [icaremedicalgroup.us](http://icaremedicalgroup.us)

## TRAINING PROGRAMS

### Educational Workshops

We now hold regular educational workshops for therapists and service providers in locations around USA.

We can provide training with our product specialists on everything from patient transfers to the psychology of care in the home.

### Product Training

Icare can provide one on one training at your showroom or facility either in person, or remotely online.

We can supply training videos and manuals, on request.





## THERAPY KITS

We have demonstration kits and portfolio's available to assist Occupational Therapist prescribe equipment.

Easy information packs to take into homes when accessing clients.



## SPARE PARTS

We carry spare parts for all our products should you need replacements.

If you are a service technician, and unsure what parts you may need, please contact us for a full parts list.



## SERVICE

In the unlikely event of a product fault or malfunction, we pride ourselves on prompt service.

If you have a product fault, or questions – please contact your dealer straight away. The dealer is responsible for servicing the beds and repairs under warranty. Patient care and comfort is very important, and that is why we do our part with quick dispatch and repairs when needed.

## FUNDING SOLUTIONS

Icare products are available through a variety of national government funding programs.

For more information on what funding options are available in your region, please contact us for your nearest dealer.

Or contact your local store to discuss your requirements.

# PARTNER BED ASSESSMENT SCRIPT FORM

## IMPORTANT!


1. Please ensure the client's bedroom layout is considered when confirming the sleeping side of the bed users (e.g. access to bathroom, carer access etc.).
2. Any two Icare bed models can be used in a partner set-up.
3. The client must sign the form to confirm they understand and confirm the bed layout before manufacturing process can begin.
4. Please complete all information on this form before submission.
5. Note: When the partner set-up (two beds) are in place, there is NO gap between the mattresses.

Client Name:  Reference:  (Must not be left blank):

Referrer Name:  Phone Number:

Email:

Equipment Supplier:  Store Contact:



Bed Size:  Bed Size:

Bed Model:  Bed Model:

Fabric Color:  Fabric Color:

Signature:

Date:



# BED PACKAGE SCRIPT FORM

Client Name:

Referrer Name:  Phone Number:

Email:

Equipment Supplier:  Store Contact:

1. Select Bed Model	2. Select Bed Size	3. Select Bed Color
<input type="checkbox"/> IC333	<input type="checkbox"/> Twin XL	<input type="checkbox"/> Onyx
<input type="checkbox"/> IC222	<input type="checkbox"/> Full XL	<input type="checkbox"/> Stone
<input type="checkbox"/> IC111	<input type="checkbox"/> Queen	<input type="checkbox"/> Smoke
<input type="checkbox"/> IC100	<input type="checkbox"/> Partner Set-up *	<input type="checkbox"/> Silver
		<input type="checkbox"/> Vinyl Ebony
		<input type="checkbox"/> Vinyl Dove

4. Select Mattress **	5. Accessories	6. Select Assistance Rails
<input type="checkbox"/> ActiveX™ IC15	<input type="checkbox"/> Headboard	<input type="checkbox"/> U-Assist
<input type="checkbox"/> ActiveX™ IC20	<input type="checkbox"/> Footboard	<input type="checkbox"/> High Side Rail
<input type="checkbox"/> ActiveX™ IC30		<input type="checkbox"/> Low Side Rail
		<input type="checkbox"/> Extendable Side Rail

7. Sheets & Protection **	8. Pillow
<input type="checkbox"/> Sheet Set - White	<input type="checkbox"/> Contour
<input type="checkbox"/> Sheet Set - Charcoal	<input type="checkbox"/> Classic
<input type="checkbox"/> Mattress Protector	<input type="checkbox"/> Curve
<input type="checkbox"/> Mattress Cover	<input type="checkbox"/> Cloud
<input type="checkbox"/> Zerotec Cover	<input type="checkbox"/> Conform

\* For Partner Set-up's it is also required to supply the 'Partner Bed Assessment Script Form'.  
 \*\* Size will be same as bed size selected.

Signature:

Date:



**ICARE MEDICAL GROUP**

**[icaremedicalgroup.com](http://icaremedicalgroup.com)**

**[us.support@icaremedicalgroup.com](mailto:us.support@icaremedicalgroup.com)**

**+1 308 351 3000**